


Cryogenic Temperature Controller

Model 22C

The Model 22C is a two-input, four-control loop cryogenic temperature controller designed for general purpose laboratory and industrial use. Each input is independent and capable of temperature measurement to <200mK with an appropriate temperature sensor. The Model 22C supports virtually any cryogenic temperature sensor produced by any manufacturer.

The four-output control loop circuits feature a primary 50W heater, a secondary heater of 25W and two non-powered outputs.

The 22C front panel incorporates a large high resolution graphics TFT type Liquid Crystal Display with an exceptionally wide viewing angle. With it's bright white LED backlight, complete instrument status can be seen at a glance, even from across the room.



Major highlights of the Model 22C:

- Two multipurpose input channels support Diode, Platinum RTD and most cryogenic NTC resistive temperature sensors. Thermocouple inputs are optional.
- Operation from 200mK to over 1500K with an appropriate sensor. Constant-Voltage, AC excitation of resistive sensors minimizes errors and extends their useful temperature range.
- Four independent control loops: Loop #1: 50-Watt, three-range; Loop #2: 25-Watt two-range; Loop #3 and #4: 10-Volt two-range.
- Large, bright and highly configurable display.
- Synchronous input filter improves control accuracy and stability in cryocooler based systems.
- Two large 10-Ampere dry-contact relays.
- Data logging to internal Non-Volatile memory.
- Table mode control automatically switches the loop input sensor to allow smooth, continuous control over a wide range of temperature.
- Remote interfaces include 100/10 Ethernet, IEEE-488.2 (GPIB) and USB 2.0. LabView™ drivers available for all interfaces.



User Interface

The Model 22C's user interface consists of a large, bright TFT type Liquid Crystal Display and a full 21-key keypad. In this user-friendly interface, all features and functions of the instrument can be accessed via this simple and intuitive menu driven interface.

```

1A:Sample      2B:Rad Shield
263.000Ω      271.321K
300.000Ω -Htr Off- 200.000K -Htr Off-
A :Sample      -- 263.000Ω --
B :Rad Shield  -- 271.321K --
Relay1: ChA
Relay2: ChB
  
```

The Home screen projects four user configurable zones that allow the real-time display of all input channel, control loop and instrument status information. From this screen, accessing any of the instrument's configuration menus requires only a single key press.

```

+   ChB:Rad. Shield
4.210 K --
Sen:20 Pt100 385
Input Conf19
CalGen
Statistics
High Alarm:200.00
High Enable:No
Low Alarm: 20.000
Low Enable:No
Deadband: 0.250
Latched Enable:No
Audible Ena: No
  
```

Innovative instrument configuration menus show real-time status information so the user can *instantly* view the results of any changes made.

```

Loop1A:Loop 1
Set Pt:300.000K   A: 0.532K
Pgain: 6.0000    1-Off-MID -Htr-Off-
Igain: 60.0005   Range: MID
Dgain: 7.5000/S  PID Table index: 1
Pman: 5.0000%    Htr Load: 50
Type: Man        *Next
Input: ChA
  
```

Flexible Inputs

The Model 22C has two identical input channels each of which implements a ratiometric AC resistance bridge.

This resistance bridge is used in a constant-voltage AC mode to provide robust support for **Negative Temperature Coefficient (NTC)** sensors including Ruthenium-oxide, Carbon-Glass, Cernox™, Carbon-Ceramic, Germanium

and several others. Since they have a negative temperature coefficient, the use of a constant-voltage measurement will reduce power dissipation in the sensor as temperature decreases. By maintaining the lowest possible power level, sensor self-heating is minimized and useful temperature range is greatly increased. An additional advantage of constant-voltage excitation is that NTC resistors lose sensitivity in the upper part of their range. By auto-ranging excitation current to maintain a constant voltage, sensitivity and noise immunity in that range is also improved.

Positive Temperature Coefficient (PTC) resistor sensors including Platinum, CLTS and Rhodium-Iron RTDs use the resistance bridge in a constant-current, AC mode. Platinum RTD sensors use a built-in DIN standard calibration curve that has been extended to 14K for cryogenic use. Lower temperature use is possible with custom calibrations.

Model 22C Supported Sensors		
	Temperature Range	Example Sensors
Silicon Diode	1.4 - 500K	Cryo-con S900, S800 SI-440, 430, 410 Lakeshore DT-670, 470
Platinum RTD	14 - 1200K	Cryo-con CP-100 Cryo-con GP-100 Cryo-con XP-100 Cryo-con XP-1K
Rhodium-Iron	1.4 - 800K	Oxford PHZ 0002
Germanium	200mK - 100K	Lakeshore GR-200A
CLTS	4 to 300K	Vishay CLTS-2B
Carbon Glass	1.4 - 325K	Lakeshore CGR-1-500
Cernox™	200mK - 325K	Lakeshore, all types
Carbon-Ceramic	200mK - 300K	TMI-A1
Ruthenium Oxide	200mK - 200K	Cryo-con R400 Cryo-con R500
Thermocouple	1.4 to 1500K	All thermocouple types.

Input Specifications					
	Diode sensors	PTC resistor sensors	NTC Resistor Sensors Constant-Current mode	NTC resistor sensors Constant-Voltage mode	Thermocouple (Option)
Input Configuration	Constant-Current DC	Constant-Current AC Resistance Bridge	Constant-Current DC Resistance Bridge	Autoranging Constant-Voltage AC Resistance Bridge	DC voltage
Input Range	0.1V - 2.25V	1.0mA: 0.1 - 500Ω 100µA: 1.0K - 5.0KΩ CLTS: 300Ω	100K - 225KΩ	Minimum: 0.5Ω Maximum: 2.0MΩ	± 70mV
Accuracy: % Rdg ± % Range	0.005% ± 80µV	0.01% + 0.0005%	±1.0KΩ @100KΩ	14Ω to 30KΩ: 0.05% + 0.05% 0.5Ω to 2MΩ: 0.15% + 0.15%	0.05%
Resolution: % Range	10µV	0.0003%	0.00004%	0.0003%	0.0003%
Excitation	10µA DC	1.0mA, 100µA	10µA DC	100mV, 10mV, 1.0mV, 300uV 2.5mA to 10nA continuously variable	N/A

Silicon diode sensors are supported over their full temperature range by using the bridge in a DC, constant-current mode.

Thermocouple sensors are supported by using an optional thermocouple module. This module plugs into any of the Model 22C's input channels and is powered by the instrument to provide amplification, cold-junction compensation and connection to copper. Up to two modules can be connected to a single instrument.

The Model 22C includes built-in **sensor calibration curves** that support most industry standard temperature sensors. Additionally, eight **user calibration curves** are available for custom or calibrated sensors. Each user curve may have up to 200 entries.

For all sensor types, conversion of a sensor reading into temperature is performed by using a **Cubic Spline** interpolation algorithm. In addition to providing higher accuracy than conventional linear interpolation, the spline function eliminates discontinuities during temperature ramps or sweeps by ensuring that the first and second derivatives are continuous.

New calibration curves may be generated using the **CalGen®** feature. This provides an easy and effective method for obtaining higher accuracy temperature measurements without expensive sensor calibrations.

Input Channel Statistics: The Model 22C continuously tracks temperature history independently on each input channel and provides a statistical summary that indicates the channel's minimum, maximum, average and standard deviation. Also shown are the slope and the offset of the best-fit straight line of temperature history data.

Data logging is performed by continuously recording input temperature data to an internal 1,365 entry buffer. Data is time stamped. Non-volatile memory is used so that data will survive a power failure.

Four Control Loops

The Loop #1 heater output is a linear, low noise RFI filtered current source that can provide up to 1.0 Ampere into 50Ω or 25Ω resistive loads. Three full-scale ranges are available in decade increments down to 500mW.

Loop #2 is a two-range linear heater with that will provide 25 or 2.5-Watts into a 50Ω load.

Loop #3 and #4 are non-powered analog voltage outputs designed to control an external booster power supply. Output is zero to 10-Volts or zero to 5.0-Volts by user selection. The controller's relays may be used to turn the external supply on and off.

All control loops are completely independent and any loop may be controlled by any sensor input. Control modes are **Manual, PID, Ramp, PID Table and Ramp Table**.

The field proven **Autotune** function of the Model 22C involves the use of a specific output waveform to first develop a process model, then generate the optimum P, I and D coefficients.

PID tables are available that can be used to store optimum control parameters vs. point temperature. Each entry of a PID table contains a setpoint, a control input, PID values and a heater output range setting. When the point is changed, the controller will automatically generate new PID values, a controlling input channel and heater range.

The Model 22C will perform a **temperature ramp** function using a specified maximum ramp rate and target setpoint.

Cryostat Protection

Damage to a cryostat or critical sample is a serious problem with any cryogenic system. The Model 22C implements the most robust set of protection features in the industry.

The **Over Temperature Disconnect** feature will disable the heater if an over temperature condition exists on any selected input channel. A fail-safe mechanical relay is used to disconnect the controller's heater thereby ensuring that the user's equipment is always protected.

The **Maximum Setpoint** feature is used to prevent the user from inadvertently entering a higher point than the equipment can tolerate and a **Maximum Power Limit** will ensure that the controller can never exceed heater power output above the set limit.

Alarms and Relays

The Model 22C supports visual, remote and audible alarms. Additionally, there are two 10-Ampere dry-contact mechanical relays. Each may be independently programmed to assert or clear based on a high or low temperature condition or a detected sensor fault.


Alarms may be latched. These are asserted on an alarm condition and will remain asserted until cleared by the user.

Remote Control

Standard Remote Interfaces include 100/10 Ethernet, IEEE-488.2(GPIB) and USB 2.0.

The Model 22C connects directly to any **100/10 Ethernet Local-Area-Network (LAN)**.

The **TCP/IP and UDP data port servers** brings fast Ethernet connectivity to all common data acquisition software programs including LabView™.



Using the Ethernet **HTTP** protocol, the instrument's **embedded web server** allows the instrument to be viewed and configured from any web browser.

LabView™ drivers are supplied for the Ethernet TCP/IP, UDP, IEEE-488.2 (GPIB) and USB 2.0 interfaces.

The Model 22C's **remote command language** is **SCPI** compliant according to the IEEE-488.2 specification. SCPI establishes a common language and syntax across various types of instruments. It is easy to learn and easy to read. The SCPI command language is identical in all Cryo-con products so that your investment in system software is always protected.

Command Scripts can be used to completely configure an instrument including setting custom sensor calibration curves and PID tables. Further, scripts can query and test data.

An **Applications Program Interface (API)** package is supplied that facilitates communication with the instrument using the TCP/IP interface. It is supplied as a Microsoft Windows™ DLL that is easily linked with C, C++ or Basic programs.


Firmware updates

Full instrument firmware updates may be installed by using the Ethernet connection. Cryo-con provides firmware updates, on request, via e-mail. They are **free of charge** and generally include enhancements and new features as well as problem fixes.

Utility Software: Utility software is provided that connects any Windows based personal computer to the Model 22C via any of its remote interfaces. This software provides a graphical control panel that greatly simplifies instrument setup and configuration. Features include:

- Continuous strip-chart monitoring of all inputs and outputs.
- Downloading, uploading, viewing and editing of sensor calibration curves and PID tables and command scripts.
- Automated instrument calibration capability.


Rear Panel Connections



- **Input Connectors:** DIN-6 recepticals provide 4-wire measurement connection plus a continuous shield through the backshell.
- **Thermocouple Option:** Connects to any of the input connectors. Up to 2 modules supported.
- **Loop #1:** 50-Watt heater output. Dual Banana Plug with chassis ground lug.
- **Loop #2:** 25-Watt heater, part of a 10 pin detachable terminal block.
- **Loop #3 and #4:** non-powered outputs. Detachable terminal block.
- **Relay #1 and #2:** Dry contact relays. Detachable terminal block.
- **Ethernet:** RJ-45 with LAN activity indicators.
- **IEEE-488.2:** Standard IEEE-488.2 (GPIB) connector.
- **USB:** USB client connector.
- **AC Power:** RFI filtered Power Entry Module including fuse drawer and line voltage selector.

Ordering Information

Product	Description
Model 22C	Controller with two multi-function sensor input channels and four control loops.
	<p>Controller includes: User's Manual, Cryo-con software CD, two input connectors, heater connector, terminal block plug, detachable power cord and a certificate of calibration.</p> <p>Specify AC Line Voltage or required power cord when ordering (may be changed in the field):</p> <ul style="list-style-type: none"> -100 Configured for 90 - 100VAC with detachable USA power cord. -120 Configured for 110 - 120VAC with detachable USA power cord. -230 Configured for 220 - 230VAC with detachable universal Euro (Shuko) line cord. -240 Configured for 240VAC with detachable universal Euro (Shuko) line cord.

Option	Description
4039-004	<p>Thermocouple Input Module. Field installable. Supports all thermocouple types.</p> 

Accessories	Description
4022-016	Input connector kit consisting of two DIN-6 sensor input connectors.
4124-018	Output connector kit consisting of a dual banana plug heater connector and a 10-pin terminal block receptacle.
4034-031	Two instrument shelf rack mount kit.
4034-032	One instrument shelf rack mount kit
4022-030	One instrument 2U rack mount kit.

Specifications

User Interface

Display Type: 240x64 graphics TFT LCD with LED back-light.
Number of Inputs Displayed: Two.
Keypad: Sealed Silicon Rubber.
Temperature Display: Six significant digits, auto-ranged.
Display Update Rate: 0.5 Seconds.
Display Units: K, C, F or native sensor units.
Display Resolution: User selectable to seven significant digits.

Input Channels

There are two input channels, each of which may be independently configured for any of the supported sensor types.

Sensor Connection: 4-wire differential. DIN-6 Connector.
Thermocouple Connection: External option. Field installable.
Sensor Types: See Supported Sensor Table.
Sensor Selection: Front Panel or remote interface.
Input Configurations: See input specifications table.
Bridge Modes: Passive, Constant-Current or Constant-Voltage.
Bridge type: Ratiometric resistance bridge. Two ADCs per input channel.
AC Excitation Frequency: 3.25Hz bipolar square wave.
Voltage Excitations: 100mV, 10mV, 1.0mV and 300uV. Minimum excitation current is 10nA, maximum is 2.5mA.
DC Offset: <8nA by active cancellation.
Sample Rate: 15Hz per channel.
Digital Resolution: 24 bits.
Measurement Accuracy: See input specifications table.
Measurement Drift: 15ppm/°C. <10Ω: 30ppm/°C.
Isolation: Input channel circuits are electrically isolated from all other internal circuitry but not from each other.
Measurement Filter: 0.5, 1, 2, 4, 8, 16, 32 and 64 Seconds.
Calibration Curves: Built-in curves for industry standard sensors plus eight user curves with up to 200 entries each. Interpolation is performed using a Cubic Spline.
CalGen®: Calibration curve generator fits any Diode or resistor sensor curve at 1, 2 or 3 user specified temperature points.

Data Logging

Data logging is performed to an internal 1,365 entry circular buffer and is time-stamped with a real-time clock. Buffer memory is non-volatile and will retain valid data without AC power. All four input channel temperatures are recorded.

User Setups

Four User Setups are available that save and restore the complete configuration of the instrument.

Control Outputs

Number of Independent Control Loops: Four.
Control Input: Any sensor input.
Loop Update Rate: 15Hz per loop.
Isolation: Control loop circuitry is referenced to chassis ground.
Control Type: PID table, Enhanced PID, Ramp or Manual.
Autotune: Minimum bandwidth PID loop design.
PID Tables: Six user PID tables available for storage of setpoint and heater range vs. PID and heater range. 16 entries/table.
Set-point Accuracy: Six+ significant digits.
Fault Monitors: Control loops are disconnected upon detection of a control sensor fault or excessive internal temperature.
Over Temperature Disconnect: Heater may be relay disconnected from user equipment when a specified temperature is exceeded on any selected input.

Loop #1 Primary Heater Output

Short circuit protected linear current source.
Ranges: Three output ranges of 1.0A, 333mA and 100mA full-scale, correspond to 50W, 5.0W and 0.5W into a 50Ω load.
Load Resistance: Selectable at 25Ω or 50Ω.
Minimum Load: 10Ω in 25W setting, 40Ω in 50W setting.
Resolution: 1.0ppm of full-scale power (20 bits).
Readbacks: Heater output power, Heatsink temperature.
Connection: Dual banana plug.

Loop #2 Heater Output

Short circuit protected linear current source.
Ranges: Two output ranges of 710mA and 224mA full-scale, which correspond to 25W and 2.5W into a 50Ω load.
Load Resistance: 50Ω for maximum output
Resolution: 1.0ppm of full-scale power (20 bits).
Readbacks: Heater output power, Heatsink temperature.
Connection: Detachable terminal block.

Loop #3 and #4 Control Outputs

Analog voltage outputs that are intended to drive an external booster supply.

Output: zero to 10 / 5.0-Volts. Output impedance: ~2,000Ω.
Resolution: 1.0ppm of full-scale power (20 bits).
Connection: Detachable terminal block.

Status Outputs

Audible and Visual Alarms: Independent audible remote and visual alarms.
Relays: Two dry-contact relays. N.O. contacts available. Ratings are 125VAC @ 10A. Maximum switching power: 150W.
Status reported via Remote Interface: Sensor fault, Heater over temperature fault.

Remote Interfaces

Maximum reading rate for all interfaces is >40 rdg/s.
Ethernet: Connects to any 10/100 Ethernet Local Area Network. Electrically isolated. **TCP/IP** and **UDP** servers provide remote control by using an ASCII command language. **HTTP** provides built-in web server. **SMTP** sends e-mail based on alarm conditions.
IEEE-488.2 (GPIB): Full IEEE-488.2 compliant.
USB 2.0: Serial port emulation. Baud Rates to 115,200.
Programming Language: IEEE-488.2 SCPI compatible.
LabVIEW™ drivers available for all interfaces.

Firmware

Instrument firmware can be updated in the field via the Ethernet connection. Firmware updates are available via the Internet free of charge.

General

Ambient Temperature: 25°C ± 5°C for specified accuracy.
Mechanical: 8.5"W x 3.5"H x 12"D. One half-width 2U rack. Instrument bail standard, rack mount kit optional.
Weight: 9 Lbs.
Power Requirement: 100, 120, 230 or 240VAC +5% -10%. 50 or 60Hz, 150VA. Field selectable.
Conformity: European CE certified.
Calibration: NIST traceable.

Contact Information

Cryogenic Control Systems, Inc.
PO Box 7012
Rancho Santa Fe, CA 92067
Tel: (858) 756-3900 Fax: (858) 759-3515
E-mail: sales@cryocon.com Web: www.cryocon.com

CalGen® and Cryo-Con® are registered trademarks of Cryogenic Control Systems, Inc. All other product and company names are trademarks or trade names of their respective companies.

©Cryogenic Control Systems, Inc. M22Cds 12/13
Specifications subject to change without notice